

Tellus 5 - Using JabblaIR with The Grid 2 and Grid 3 for Environmental Control


This guide provides information on setting up environmental control on the Tellus 5 for Grid 2 and Grid 3.

The environmental control unit contained in the Tellus 5 is called JabblaIR. JabblaIR works in a similar way to other environmental control modules in that it allows you to record and send infrared (IR) signals to control appliances such as televisions, DVD players and infra-red remote controlled toys. Recorded infrared signals can then be transmitted from the Tellus 5.

On the Tellus 5 the IR receiver is located on the front of the device in the top right hand corner. The IR receiver is where you will point the remote control that you are trying to record.



The Tellus 5 IR transmitter is located on the back of the device. Like all IR remote controls, this transmitter must have a direct line of sight to the device (e.g. TV) that you want to



When recording remote controls it is highly recommended to have brand new, good quality batteries installed in the remote control.

Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353



Recording Remote Control Commands in Jabbla IR


To record remote controls on the Tellus 5 for use with Grid software, a small program called JabblaIR Tool is used. The tool can be found on the desktop of the Tellus 5 or alternatively can be downloaded from Zyteq's website.

To record a remote:


Open Jabbla IR Tool.
 Select Run if asked by Windows.




2. Select Add New (green + sign).




- 3. In the **Name** box type the name of the remote code being recorded.
 - Be sure to be clear with the name so that you are able to easily identify which command has been recorded. E.g. TV Power. It is also recommended that you also complete the IR record on the final page of the instructions.
- 4. Once you have finished editing the name of the remote code select **Record.**
- 5. Point the remote control at the IR Reciever (on the front of the device) from approximately 15-20cm away and press the button on the remote that you want to record. Once you have sucessfully recorded the remote the **Status** display will say **OK** and turn green.
- 6. Press **Test** to check that the recording worked (make sure that the IR transmitter is facing the TV/DVD etc).










Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353




To use a specific remote code that has been recorded in JabblaIR you need to be able to create some different arguments. An argument is a small amount of text that tells the JabblaIR tool which remote code to activate and how to activate it. JabblaIR can use a single recorded remote button or a combination of buttons (eg. Chanel 24 is accessed via the number 2 and the number 4). This kind of remote activation where multiple buttons are automatically used is called a macro. Some macro commands will also require the use of a wait function for devices that require a gap between button presses.

An argument can contain the following types of remote activation:



Below are some examples of types of arguments that can be created. In these examples the recorded remote buttons were as pictured.



Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353



There are four types of arguments that may be used.

- 1. Simple argument a single remote code is activated
 - In this case the argument is simply the number of the of the remote code in Jabbla IR.
 - o Eg. The argument for Chanel Up to be activated once by itself is simply 4.
- 2. Repeating argument a single remote code is activated a number of times.
 - o In this case a single remote button is pressed repeatedly.
 - Eg. The argument for Volume Up pressed 5 times is 2[space]5 Where the
 first number is the remote code to be activated and the second number is
 the times it is repeated.
- 3. Simple macro a number of different buttons are pressed automatically.
 - In this case different buttons are pressed automatically, one after the other. This can be useful for changing to specific TV channels.
 - Eg. The argument for channel 12 would be 6 + 7 where the sixth IR command (Number 1 on remote) is pressed and then the seventh IR command (Number 2 on remote) is pressed immediately after.
- 4. Complex macro multiple commands are pressed with wait times.
 - In this case multiple different buttons are pressed (some repeatedly) with additional wait times between remote presses.
 - Eg. The argument for changing the channel to 12 plus plus some wait time and then turning the volume up 5 times would be 6 + 7 + s1000 + 2[space]5



Using Remote Control Commands with Grid 2


 In editing mode select the cell that is going to have the environmental control command. Select Edit cell.



2. Select add a command.



3. Search for **Run Program**. Select the command and then press **OK**.



4. Select Browse...




ZYTEQ pty ltd PO Box 190 South Melbourne Victoria 3205


Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353



Navigate to the **Desktop** and select **JabblaIR Tool.** Select **Open.**



- 6. In the text field scroll using the right arrow on the keyboard to the end of the line.
- 7. At the end of the line type a single SPACE and then type the argument.
 Eg. From the previous example where the argument to change the TV channel to 12 would be 6 + 7 the line of text would look like:
 - C:\Users\Desktop\JabblaIRTool.exe 6 + 7



8. Select OK.


Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353




Using Remote Control Commands with Grid 3:

1. In editing mode select the cell that is going to have the environmental control


command. Select Add Command



2. Search for **Start Program**. Select the command and then select **Next**.



3. Select Browse.




ZYTEQ pty ltd PO Box 190 South Melbourne Victoria 3205


Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353



4. Navigate to the **Desktop** and select **JabblaIR Tool.** Select **Open.**



5. Select Jabbia IR Tool and then select OK



6. In the **Arguments** box type the argument. Eg. From the previous example where the argument to change the TV channel to 12 would be **6 + 7** the command would look like:



Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353







Recording List

Number	Name	Location	Device	Function	Notes
eg. 1	TV Power	Bedroom	Panasoníc TV	Power On/OFF	

Ph: 03 9696 2944 Fax: 03 9696 1755 Freecall: 1800 818 353