

# BIGtrack 2

Product Number: 12000006

MSRP: ~~\$100.00~~

**\$70.00**

[See More Info](#)

Quantity

1

Add to Cart

*Free standard shipping on all orders & most orders ship within 1-business day.*

Make sure you get the correct product by providing your purchaser the product number and product name. Click to print or export a PDF to provide your purchaser with this information.



Roll over image to zoom in



## About This Product

### Overview

BIGtrack 2 is a USB mouse that includes a large 3-in / 7.62-cm trackball. The sturdy base of BIGtrack 2 allows the user to easily move and position the mouse cursor on the computer screen while using fine or gross motor movements of the hand, arm, or foot. Once the mouse cursor is positioned, two over-sized buttons are available for left and right clicking.

#### **BIGtrack 2 Includes:**

- Large 3-in / 7.62-cm trackball
- Updated design
- Two switch jacks for left and right click
- Color differentiated right and left click buttons
- USB connection (USB Type A receptacle required)

#### **Drag Lock**

The drag lock feature enables a user to move objects on the screen without having to hold the mouse click button down. To activate the drag lock feature, position the mouse cursor over the item you would like to move. Once the cursor is positioned, double right click to lock onto the item. Once locked, the user can move the item with the trackball. Once a desired location is reached, activate either the right or left click button one time to release the object.


#### **Switch Access**

BIGtrack 2 includes two 1/8-in / 3.5-mm mono (TS) switch jacks. Plugging in a switch on the left side of BIGtrack 2 will enable the left click via the switch. Plugging in a switch on the right

side of BIGtrack 2 will enable the right click via the switch. Switch are not included and sold separate.

## Compatibility

BIGtrack 2 is compatible with most laptop, desktop, and tablet computers that can use a standard mouse. One USB Type A receptacle is required.

 Make sure you get the correct product by providing your purchaser the product number and product name. Click to print or export a PDF to provide your purchaser with this information.

---

## ableCARE, ableSAFE, and ableSTRONG



Our ableCARE program removes stress and barriers in implementing AbleNet solutions. The ableCARE Product Success team is an extension of your school or clinic team and is standing by ready to provide guidance and knowledge so you and your team build confidence in implementing the AbleNet solutions you've invested in.



Certified third-party laboratories test ableSAFE products to ensure they meet or exceed current global safety standards. When you see the ableSAFE badge, you know you are using a device engineered with safety in mind, made from high-quality materials, and safe for all users.



ableSTRONG products are designed to be dropped, encounter moisture and dirt, and keep performing so users can achieve best outcomes. Before being approved as ableSTRONG, AbleNet products go through a series of durability tests to ensure they can stand up to a variety of environmental factors.

[Learn More](#)

---

## FAQs and Downloads

[Visit our knowledge base](#) for product questions, getting started guides, instructions for use, and other downloads.

---

## Technical Specifications

- Trackball size: 3-in (7.6-cm)
- Connection type: USB and PS2 Adapter

## Regulatory


This product is compliant with Medical Device Regulation 2017/45. [Download the declaration of conformity.](#)

## Recently Viewed



Buddy Button

**\$52.50** MSRP: ~~\$75.00~~




LITTLEmack

**\$124.00** MSRP: ~~\$155.00~~



BIGtrack 2

**\$70.00** MSRP: ~~\$100.00~~



BIGmack

**\$124.00** MSRP: ~~\$155.00~~



BIG Step-by-Step

**\$176.00** MSRP: ~~\$220.00~~

[Manage Website Data Collection Preferences](#)

Welcome to AbleNet! 🙌

Have a question? The AbleNet team is here to help!

[Chat with the team now](#)